



[http:// www.femtopath.com](http://www.femtopath.com)

*FemtoPath*

EGFR BRAF AKT APT NOTCH1 NRAS  
JAK2 KIT KRAS CTNBT Her2 PLS  
PTEN

*Make molecular diagnostic easier*

# CYP2C19 Genotyping HRM Kit

Catalog Number : CYP-001-0020

20 tests/kit

For Research Use Only

Website : [www.femtopath.com](http://www.femtopath.com)

E-mail : [info@hongjing.com.tw](mailto:info@hongjing.com.tw)

Address : 5F., No.172, Sec. 1, Zhongshan Rd., Yonghe Dist., New Taipei City 234, Taiwan (R.O.C)

Tel : +886 2 32338585

Fax : +886 2 32338686



## Background Information

Cytochrome P450 2C19 (abbreviated CYP2C19) gene produces an enzyme which is involved in protein processing and transport in liver cells. CYP2C19 enzyme plays a role in the processing or metabolizing of at least 10 percent of commonly prescribed drugs, including a drug called clopidogrel (also known as Plavix). Clopidogrel is an antiplatelet drug that CYP2C19 enzyme activates its active form to inhibit receptor protein on the surface of platelets. CYP2C19\*2 and CYP2C19\*3 are most common CYP2C19 gene polymorphisms associated with clopidogrel resistance result in the production of a nonfunctional CYP2C19 enzyme that is unable to activate clopidogrel. In contrast, CYP2C19\*17 polymorphism enhances antiplatelet activity of clopidogrel. CYP2C19 also play a important role in the metabolism of Omeprazole. Omeprazole is Proton pump inhibitors (PPIs) that are used in the treatment of stomach-acid related disorders. The CYP2C19\*17 variant allele which is associated with increased enzyme activity may enhance PPI clearance.

FemtoPath CYP2C19 Genotyping HRM Kit is a highly specific and sensitive HRM-based test able to detect genetic variants in CYP2C19. Femtopath CYP2C19 Genotyping HRM Kit enables to identify CYP2C19\*2, CYP2C19\*3 and CYP2C19\*17 variants. CYP2C19\*2, CYP2C19\*3, and CYP2C19\*17, could explain most of the phenotypes related to the voriconazole metabolism and some of its pharmacokinetic singularities.

## Intended Use

Identification of CYP2C19 genetic variants including CYP2C19\*2, CYP2C19\*3 and CYP2C19\*17.

## Kit Components

Vial Number	Name of Component	Description	Volume (µl)
1	CYP2C19*2 Primer Mix	Position 681 of the gene (G/A)	80
2	CYP2C19*3 Primer Mix	Position 636 of the gene (G/A)	80
3	CYP2C19*17 Primer Mix	Position -806 of the gene (C/T)	80
4	2X HRM Master Mix	qPCR reaction polymerase pre-mixture	200 x 3
5	DDH2O	Sterile H2O	1000



## Additional Equipment and Reagents Not Supplied

DNA extraction reagents and equipments  
Sterile, nuclease-free PCR tubes  
Pipettes and appropriate pipette tips  
qPCR instruments

## Storage

Store at -20°C. Stable for one year from the date of shipment. Non-hazardous. No MSDS required.

## DNA Preparation

Human genomic DNA must be extracted from formalin-fixed paraffin-embedded (FFPET), fine needle biopsy or pleural effusion specimens. We recommend using Qiagen DNA extraction kit (QIAamp DNA FFPE Tissue Kit, cat No. 56404) for genomic DNA extractions. The OD value of genomic DNA extractions should be measured using the spectrophotometer. Extracted genomic DNA specimens may be stored at -20 °C for long term storage.


## FemtoPath CYP2C19 Genotyping HRM Test Preparation

1. Thaw and centrifuge all vials in the kit.
2. Keep 2x HRM Master Mix in ice and dark to prevent enzyme degradation
3. Prepare separately qPCR Reaction Mixture by adding 10 µl 2x HRM Master Mix, 10 µl CYP2C19 Primer Mixes and 5µl Sterile H<sub>2</sub>O with a total of 19 µl mixture per reaction.
4. Add 1 µl (20~100 ng) DNA specimen into the PCR reaction mixture.

Components	Volume (µl)
DNA specimen (20~100 ng)	1
CYP2C19 Primer Mixes	4
2x HRM Master Mix	10
DDH <sub>2</sub> O	5
Total Volume	20

## QuantStudio™ 3 or 5 Real-Time PCR Cycling Condition

Centrifuge all PCR mixture and perform qPCR using the cycling conditions below.




## Assign Targets and Samples

Quick Setup **Advanced Setup**


Targets		Add		Action			
	Name	Reporter	Quencher	Comments	Task	Quantity	
<input type="checkbox"/>	Target 1	Eva dye	None		▼		✕

## Data Analysis


CYP2C19\*2



CYP2C19\*3



CYP2C19\*17





## Warnings and Precautions

1. Please read the instruction carefully before use.
2. Experiments should be performed under proper sterile condition with aseptic techniques.
3. Thawed thoroughly, mix the components by inverting and centrifuge briefly before use.
4. Use sterile pipette tips and dedicated pipettes to prevent exogenous DNA contamination .
5. Strictly not recommended to add other sources of reagents.

## References

1. Kim, K. A., Park, P. W., Hong, S. J. & Park, J. Y. The effect of CYP2C19 polymorphism on the pharmacokinetics and pharmacodynamics of clopidogrel: a possible mechanism for clopidogrel resistance. *Clin Pharmacol Ther* 84, 236-242, doi:10.1038/clpt.2008.20 (2008).
2. Mega, J. L. et al. Cytochrome p-450 polymorphisms and response to clopidogrel. *N Engl J Med* 360, 354-362, doi:10.1056/NEJMoa0809171 (2009).
3. Collet, J. P. et al. Cytochrome P450 2C19 polymorphism in young patients treated with clopidogrel after myocardial infarction: a cohort study. *Lancet* 373, 309-317, doi:10.1016/S0140-6736(08)61845-0 (2009).
4. Holmes, M. V., Perel, P., Shah, T., Hingorani, A. D. & Casas, J. P. CYP2C19 genotype, clopidogrel metabolism, platelet function, and cardiovascular events: a systematic review and meta-analysis. *JAMA* 306, 2704-2714, doi:10.1001/jama.2011.1880 (2011).
5. Chen, K. et al. Impact of the CYP2C19 Gene Polymorphism on Clopidogrel Personalized Drug Regimen and the Clinical Outcomes. *Clin Lab* 62, 1773-1780, doi:10.7754/Clin.Lab.2016.160216 (2016).
6. Charfi, R. et al. Response to clopidogrel and of the cytochrome CYP2C19 gene polymorphism. *Tunis Med* 96, 209-218 (2018).
7. Baldwin, R. M. et al. Increased omeprazole metabolism in carriers of the CYP2C19\*17 allele; a pharmacokinetic study in healthy volunteers. *Br J Clin Pharmacol* 65, 767-774, doi:10.1111/j.1365-2125.2008.03104.x (2008).
8. Dean, L. in *Medical Genetics Summaries* (eds V. M. Pratt et al.) (2012).
9. Kamiya, C. et al. Effect of co-administered inducer or inhibitor on omeprazole pharmacokinetics based on CYP2C19 genotype. *J Pharmacol Sci* 139, 361-366, doi:10.1016/j.jphs.2019.03.001 (2019).

## Contact Information

HongJing Biotech Inc.,

Address : 5F., No.172, Sec. 1, Zhongshan Rd., Yonghe Dist., New Taipei City 234, Taiwan (R.O.C)

Phone : +886-2-3233-8585

Fax : +886-2-3233-8686

**Technical Support E-mail : [info@hongjing.com.tw](mailto:info@hongjing.com.tw)**